

LABA - Configuration

Configuration

- [Lab Test Templates](#)
- [Priority setup](#)
- [Parameters](#)
- [Treatments](#)

Lab Test Templates

Lab test templates are predefined set of tests templates setup

Tests are grouped for various industrials

User can setup the group that they want to view

Lab



Configuration

Template Tests

Parameters

Treatments

Priority

Template Tests search

Search Edit

General

Description: Test Group:

Page 1 of 1

#	Seq No	Description	Test Group	Quick Add	Active	Short Name
1	1	MD Germination	Seed Tests	Yes	Yes	MD-G
2	2	TZ	Seed Tests		Yes	TZ
3	3	QS Germination	Seed Tests		Yes	QS-G
4	4	UV	Seed Tests		Yes	UV
5	5	Seed Count	Seed Tests		Yes	Seed
6	6	Moisture	Seed Tests		Yes	Mois
7	10	Endophyte Seed	Seed Tests		Yes	Endo
8	7	Purity	Seed Tests		Yes	Pure
9	8	FD Germination	Seed Tests		Yes	FD-G
10	9	KASP	Seed Tests		Yes	KASP
11	11	Endophyte Seedling	Seed Tests		Yes	Endo-Sdling
12	12	Endophyte Tiller	Seed Tests		Yes	Endo-Till
13	13	Endophyte QS	Seed Tests		Yes	Endo-QS
14	14	Phytosanitary Inspection	Seed Tests		Yes	PHY-INS

Test Types	Details
Search button	Search the Predefined tests
Edit	Edit the test to put recommended, Test parameters in Can change the Test Description
Test Description	Search the grid for a specific test
Test Group	Dropdown of predefined industry tests

Edit Test Template

SEED001 - MD Germination test

There are three types

Tabs defined	Details
Details	Standards for each test with defaults
Test	The test varies according to the test

×

Test Template Edit

Detail

Test

Invoice

Description:

MD Germination

Short Name:

MD-G

Test Group:

Seed Tests

Test:

SEED001

Sequence:

1

Colour:

Blue

Parameter:

Parameter one set 3

Recommended:

test

Treatments:

Treatment 1

Treatment 2

Treatment 3

Add Treatments

Active:

☒

Show on Quick Add:

☒

Save

Save and Close


Close

Test Type	Details
SEED001	MD Germination test
Test Description	Description of test - this can be edited
Active	Test is active or inactive
Test Type	Predefined tests - This cannot be edited, but shows as to which name is given to test
Test Group	Predefined industry group of tests, This cannot be edited, but shows which industry the test belongs to
Sequence	The sequence of test to be shown on the quick add test

Show on Quick Add	Indicate whether this test will show as a list item on the quick add test
Recommended	Add in text to show when the test is added to a Lab test
Test Parameters	Add in text to show when the test is added to a Lab test
Treatments	Add in text to show when the test is added to a Lab test
Test results	Sample of what test results view will look like when added to the Lab Test

Add Treatments button

This adds in the Treatments description into the Treatment field



The test tab

Each test will have its own defined layout

×

Test Template Edit

Detail

Test

Invoice

Wt Per 100

A	B	C	D	E	F	G	H	Total	1000 Seed Weight
0	0	0	0	0	0	0	0	0	undefined

Save

Save and Close

Close

Invoice Tab

Define the Service or custom item to invoice for test when invoicing

Test Template Edit

Detail

Test

Invoice

Item:

BAT/BOX

Item Description:

Battery Box

Quantity:

1

Priority setup

Setting up the Priority types

Lab



Configuration

Template Tests

Parameters

Treatments

Priority

Priority Search

Search

Add

Edit

Delete

Print Grid

Code:

Description:

Active: ☒

Priority

Page 1 of 1

#	Seq No	Icon	Code	Description	Active	Is Default
1	2		PriorityA	Top Priority	Yes	
2	3		NORM	Normal Priority	Yes	Yes
3	3		PC	Priority C %, ", '	Yes	

Add / Edit

Priority Add / Edit

Code:

NORM

Description:

Normal Priority

Icon:

Circle Green

Sequence:

3

Active:

☒

Default:

☒

Save and Close

Close

Priority Edit	Details
Code	Enter a Unique code for the Priority
Description	Enter a description for Priority
Icon	Select an icon from dropdown list
Active	Make this item active
Default	Only one item can be a default Used when adding a new Lab test

Parameters

Setup default Parameters

Search for Parameter

Lab



Configuration

Template Tests

Parameters

Treatments

Priority

Parameter search

Search

Add

Edit

Delete

Print Grid

Parameter Import Export

Code:

Description:

Active: ☒

Parameters

Page 1 of 1

#	Code	Description	Active
1	PAR01	Parameter one set 3	Yes
2	PAR02	Parameter two set	Yes
3	PAR03	Parameter one set 4	Yes
4	PAR04	Parameter one set %, ", '	Yes
5	Temperature 20	Temperature 20 degrees	Yes

Add / Edit Parameter

Parameter Add / Edit

Code:

Description:

Active: ☒

Save and Close

Close

Treatments

Setup treatments

Search Treatments

Lab



Configuration

Template Tests

Parameters

Treatments

Priority

Treatment search

Search Add Edit Delete Print Grid Treatments Import / Export

Code: Description:

Active: ☒

Treatment

Page 1 of 1

#	Code	Description	Active
1	KNO3	KNO3	Yes
2	PRECHILLED	Test Prechilled 5°C	Yes

Add / Edit Treatments

Treatment Add / Edit

Code: PRECHILLED

Description: Test Prechilled 5°C

Active: ☒

Save and Close Close